

PROFESSIONAL ROOM STERILIZERS


STERYLIS®

ULTRA



PROFESSIONAL **ROOM** STERILIZERS

Sterilization by ozone and UV-C radiation with particle filtration in one unit.



TECHNOLOGY



ULTRA ROOM STERILIZERS

These are innovative room disinfection devices with an additional air filtration function, which use UV-C sources with the wavelength of the most effective in the fight against microorganisms, and efficient ozone generators that enable a quick and effective disinfection process at even hard-to-reach places. They are available in several versions differing in the power of the UV-C source used, the efficiency of the ozone generator and the airflow efficiency. The series include the models: ULTRA-220, ULTRA-330, ULTRA-440, ULTRA-550.



UV-C RADIATION

The high-performance lamp system between the filters emits UV-C radiation with wavelength of 253.7 nm. UV-C radiation causes irreversible damage to the DNA and RNA of microorganisms (bacteria, viruses, molds, fungi) leading to their neutralization [13] [14]. This means that the air after passing through the channel contains a significantly reduced amount of active microorganisms. Thanks to the special closed disinfection channel construction, the emitted high-energy UV-C radiation does not spread outside the interior of the unit, which allows for safe operation of the sterilizer in this mode in rooms where people are present.



OZONE GENERATOR


In the intensive disinfection mode, the device produces ozone, which has strong oxidizing properties – including the destruction of structures of living organisms [6] [7] [8] [9] [11] [12]. Ozone also has the ability to neutralize all kinds of odors. Thanks to the high performance of the ozone generators used in the unit, the speed and efficiency of this method of decontamination with STERYLIS sterilizers is extremely high. The total sterilization of the rooms is usually completed in about 4 hours, and the whole process is controlled by a controller measuring the ozone concentration generated in the room in real time. Intelligent control of the ozone process ensures not only the highest sterilization efficiency but also the safety of the entire process. The safe concentration of ozone is possible thanks to the function of its destruction after the sterilization process. Thanks to the gaseous form of the ozone generated, not only the air in the room is completely sterilized, but also all the objects in the room that ozone can reach during the decontamination process.



DOUBLE FILTRATION

STERYLIS sterilizers have 2-stage air filtration (pre-filter and exhaust filter). High-quality filter media with electrostatic properties are capable of trapping even microscopic dirt particles. The use of an antibacterial system based on an odorless, non-toxic and invisible to the human eye coating applied to the filter medium allows to remove harmful allergens and bacteria from the air. This technology makes the STERYLIS sterilizers, in addition to their disinfection and sterilization functions, extremely effective also in cleaning air from particulate matter.

INTUITIVE AND CONVENIENT



1. **LIGHT AND SOUND SIGNALLING**
indicates that the unit is in the room
sterilization cycle with ozone.
2. **ERGONOMIC HANDLES**
3. **CONTROL PANEL**
4. **MAIN SWITCH**
5. **CARTRIDGE FILTERS**
easy to replace
6. **CONTROLLER WITH OZONE SENSOR**
manages the operation of the device and
analyzes the concentration of ozone in the
room
7. **BACKLIT PANEL**
gently pulsating blue light
8. **CASTORS**

STERYLIS®

PROFESSIONAL ROOM STERILIZERS



SAFETY

Signal and warning lights and acoustic signals generated by the device inform the user about the sterilization mode activated and the safe or excessive level of ozone concentration in the room. Whenever the device detects that the concentration limit has been exceeded, it adjusts its own performance while ensuring maximum safety.



CONVENIENCE

Intelligent air sterilizers provide maximum comfort. A simple and intuitive control panel makes it easy to select the operating mode. The new sterilization mode controller allows you to choose between automatic operation and operation with predefined settings for sterilization time and ozone concentration, it is also possible to program and select own operation mode. The displays inform the user about the ozone concentration and process parameters. They indicate the need to replace filters, the need to replace UV-C lamps if their service life is exceeded or if they burn out, and damage to other components, e.g. one of the ozone generators, if such damage occurs. This is by far the highest level of self-diagnosis in such units.



SIMPLE OPERATION

The controller the unit is equipped with is intuitive, functional and very simple to use. The operating modes are changed by means of a mode selection buttons. Signal lights visible from a distance allow you to conveniently assess from a distance which operating mode the device is currently in; safety messages are signaled in the same way.


The panel has displays indicating the current ozone concentration, time to process completion, in case of delayed activation, operating modes or selected program, maintenance alerts for the device, including signaling the need to replace filters. Along with acoustic signaling, the flashing signal lamp visible from a distance provide safety in sterilization mode. These are just some of its functions.




SILENT OPERATION

By changing the operating mode to "UV-C disinfection SILENT MODE", the unit switches to low fan speed operation and starts quiet operation in disinfection mode. From now on, the user's peaceful sleep with the simultaneous operation of the sterilizer is not at risk.


HOW THE UV-C TECHNOLOGY WORKS FOR PATHOGENS




- 1 Bacteria, viruses, and other pathogens are transferred to the sterilizer.




- 2 Pre-filtration of particulates.



- 3 Pathogens in the air are exposed to UV-C radiation.



- 4 The DNA / RNA structure of pathogens is damaged, which prevents its reproduction.



- 5 The cleaned air is returned to the room through a second outlet filter.



WHAT A FULL CYCLE OF OZONE STERILIZATION LOOKS LIKE

1

The molecular oxygen contained in the air is supplied to an ozone generator located inside the sterilizer.



2

Crown discharges separate oxygen molecules.




3

Single oxygen atoms bond with oxygen molecules to form ozone molecules (O_3).




4

Once the concentration of ozone as a gas has been achieved, it reaches all corners of the room, neutralizing pathogens.



5

Once the sterilization process is complete, the unit enters into UV ozone destruction mode, speeding up the time needed to achieve a concentration of O_3 safe for the user.




6

Once the safe ozone concentration is reached, the sterilization cycle is completed and the sterilizer enters the standby mode.



7

The room is sterile, free of pathogens.



UV-C RADIATION DOSE

GENERATED BY THE EQUIPMENT

DURING DISINFECTION

The degree of deactivation of pathogens with ultraviolet radiation directly depends on the dose of this UV-C radiation. The UV-C dose is the product of the irradiance [I] and the exposure time [t].




Therefore: DOSE = I x t

The UV-C radiation dose is expressed in joules per square meter [J/m²].

The average dose of UV-C radiation generated at one passage of air through the full volume of the disinfection chamber



| STERYLIS unit model | | ULTRA-220 | ULTRA-330 | ULTRA-440 | ULTRA-550 |
|--|-----------------------|-----------|-----------|-----------|-----------|
| Operation mode UV-C disinfection (silent) | Average dose [J/m²] | 144 | 221 | 105 | 104 |
| | Efficiency [m³/h] | 170 | 185 | 580 | 680 |
| Operation mode Disinfection UV-C | Average dose [J/m²] | 106 | 141 | 90 | 79 |
| | Efficiency [m³/h] | 230 | 290 | 680 | 900 |



UV-C RADIATION DOSE

REQUIRED TO NEUTRALIZE THE PATHOGEN

The UV-C radiation dose required for a 99.9% pathogen reduction is shown in the table. The exemplary data given are collected from scientific publications and research on ultraviolet light sterilization technology from around the world. [1] [2] [3] [4] [5] [8] [10] [12]

| BACTERIA | [J/m²] |
|--|----------|
| Agrobacterium Lumefaciens | 85 |
| Acinetobacter | 14 |
| Bacillus Anthracis | 87 |
| Bacillus Anthracis Spores | 462 |
| BaciHus Megatherium Sp. (Veg) | 25 |
| Bacillus Megatherium Sp. (Spores) | 52 |
| Bacillus Paratyphosus | 61 |
| Bacillus Subtilis | 110 |
| Bacillus Subtilis Spores | 220 |
| Bordetella pertussis | 63 |
| Clostridium Tetani | 231 |
| Clostridium Botulinum | 112 |
| Corynebacterium Diphtheriae | 65 |
| Dysentery Bacilli | 42 |
| Eberthella Typhosa | 41 |
| Enterobacter cloacae | 64 |
| Enterococcus | 28 |
| Escherichia Coll | 86 |
| Haemophilus influenzae | 19 |
| Haemophilus parainfluenzae | 77 |
| Klebsiella pneumoniae | 52 |
| Legionella pneumophila | 11 |
| Legionella Dumoffill | 55 |
| Legionella Gormanil | 49 |
| Legionella Micdadei | 31 |
| Legionella Longbeachae | 29 |
| Legionella Pneumophfla | 27 |
| Leptospiracanicola - Infectious Jaundice | 60 |
| Leptospira Interrogans | 80 |
| Micrococcus Candidus | 123 |
| Micrococcus Sphaeroides | 154 |
| Mycobacterlum Tuberculosis | 100 |
| Neisseria Catarrhalis | 85 |
| Phytomonas Tumefaciens | 105 |
| Proteus Vulgaris | 39 |
| Pseudomonas Aeruginosa | 66 |
| Pseudomonas Fluorescens | 76 |
| Rhodospirillum Rubrum | 61 |
| Salmonella Enteritidis | 100 |
| Salmonella Paratyphi | 152 |
| Salmonella Species | 70 |
| Salmonella Typhimurium | 105 |
| Salmonella Typhosa | 264 |
| Salmonella | 61.6 |
| Sarcina Lutea | 42 |
| Serratia Marcescens | 34 |
| Shigella Dysenteriae-Dysentery | 34 |
| Shigella Flexneri-Dysentery | 70 |
| Shigella Paradyntenteriae | 85 |
| Shigella Sonnei | 66 |
| Spirillum Rubrum | 61.6 |
| Staphylococcus Albus | 57.2 |
| Staphylococcus Aureus | 66 |
| Staphylococcus Epidermidis | 58 |
| Streptococcus Faecaila | 100 |
| Streptococcus Hemolyticus | 55 |
| Streptococcus Lactis | 88 |
| Streptococcus Pyrogenes | 42 |
| Streptococcus Salivarius | 42 |
| Streptococcus Viridans | 38 |
| Typhoid Fever | 41 |
| Vibrio Comma (Cholera) | 65 |
| Vibrio Cholerae | 65 |


| MOLDS | [J/m²] |
|-------------------------|----------|
| Aspergillus Amstelodami | 770 |
| Aspergillus Flavus | 990 |
| Aspergillus Glaucus | 880 |
| Mucor Mucedo | 770 |
| Mucor Racemosus (A & B) | 352 |
| Oospora Lactis | 110 |
| Penicillium Chrysogenum | 560 |
| Penicillium Digitatum | 880 |
| Penicillium Expansum | 220 |
| Penicillium Roqueforti | 264 |

| PROTOZOANS | [J/m²] |
|----------------------------|----------|
| Chlorella Volgaris (atgae) | 220 |
| E. Hystolytica | 840 |
| Nematode Eggs | 400 |

| VIRUSES | [J/m²] |
|-----------------------------|----------|
| Adenovirus | 45 |
| Bacteriophage (E.Coli) | 66 |
| Coronavirus (SARS) | 18 |
| SARS Coronavirus CoV-P9 | 40 |
| Murine Coronavirus (MHV) | 103 |
| SARS Coronavirus (Hanoi) | 134 |
| SARS Coronavirus (Urbani) | 241 |
| Coxsackievirus | 63 |
| Infectious Hepatitis | 80 |
| Influenza | 34 |
| Measles virus | 22 |
| Mumps virus | 30 |
| Norwalk virus | 198 |
| Parainfluenza virus | 21 |
| Parvovirus B19 | 25 |
| Poliovirus | 210 |
| Reovirus | 158 |
| Rhinovirus | 162 |
| Rotavirus | 240 |
| RSV | 25 |
| Rubella virus | 622 |
| VZV (Varicella surrogate k) | 18 |
| Variola | 240 |

| FUNGI | [J/m²] |
|---------------------------------|----------|
| Aspergillus spores | 258 |
| Baker's Yeast | 88 |
| Blastomyces dermatitidis spores | 140 |
| Brewer's Yeast | 66 |
| Common Yeast Cake | 132 |
| Cryptococcus neoformans spores | 138 |
| Fusarium spores | 269 |
| Mucor spores | 228 |
| Rhizopus spores | 267 |
| Saccharomyces Cereisiae | 132 |
| Saccharomyces Ellipsoideus | 132 |
| Saccharomyces Sp. | 176 |

WHAT THE OZONE PROCESS IN STERILIZED ROOMS LOOKS LIKE




Key:
 T_s - time to reach the sterilizing ozone concentration (default 5 ppm)
 T_o - supervised by the controller, required time of the sterilization process depending on the achieved ozone concentration
 T_d - supervised by the controller, supported by an ozone destructor, time of the process of returning to a safe ozone concentration

Application range
in sterilization mode (O_3)

| MODEL | RECOMMENDED STERILIZED CUBIC CAPACITY (in ozonization mode) [m³] | MAXIMUM STERILIZED CUBIC CAPACITY (in ozonization mode) [m³] |
|-----------|--|--|
| ULTRA 220 | 450 | 2,700 |
| ULTRA 330 | 850 | 3,600 |
| ULTRA 440 | 850 | 3,600 |
| ULTRA 550 | 1,050 | 5,600 |


Chart of ozone sterilization time depending
on the volume of the room to be sterilized



PLACES OF APPLICATION

PRACTICALLY IN EVERY ROOM!

STERYLIS ULTRA devices are adapted to every type of room and to different cubic capacities. Their power and performance allows them to be adjusted to operate both during the day and at night (night mode - silent). They are characterized by many advantages, such as safety, convenient usage / service / and simple operation.



Recommended
CUBIC CAPACITY
of the sterilized room

450 m³

WORK MODES:



Filtration



Disinfection
UV-C



UV-C
disinfection -
silent mode




Sterilization
O₃ + UV-C (Auto)



TECHNICAL SPECIFICATIONS:

| | |
|---|-----------------------|
| MODE – STERILIZATION O₃ + UV-C | |
| Recommended maximum cubic capacity of the sterilized room (to achieve 5 ppm O ₃) | 450 m ³ |
| Ozone generator capacity | 20,000 mg/h |
| Fan capacity Sterilization mode O ₃ +UV-C | 290 m ³ /h |
| Noise level – sterilization mode O ₃ +UV-C | 52 dB(A) |
| STANDARD MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 483 J/m ² |
| Average UV-C radiation dose | 106 J/m ² |
| Fan performance UV-C disinfection mode | 230 m ³ /h |
| Noise level – disinfection mode UV-C | 50 dB(A) |
| SILENT MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 654 J/m ² |
| Average UV-C radiation dose | 144 J/m ² |
| Fan performance UV-C disinfection mode (silent) | 170 m ³ /h |
| Noise level – disinfection mode UV-C (silent) | 44 dB(A) |
| FILTRATION MODE | |
| Fan capacity in filter mode | 230 m ³ /h |
| Noise level – filtering mode | 50 dB(A) |
| UV lamp | |
| | YES |



| | |
|------------------------------|---|
| Type lamps UV | UV-C sterilization λ = 253.7 nm |
| Electrical power of UV lamps | 220 W |
| Durability of UV lamps | 9,000 h |
| Power of the UV light source | 78 W |
| Ozone generator | YES |
| Ozone concentration sensor | YES |
| Ozone destructor | YES |
| Air filtering | 2-stage |
| Fan type | Radial |
| Fan engine | enclosed |
| Type | flowtype |
| Controller | YES |
| Automatic mode of operation | YES (control of sterilization time according to actual ozone concentration readings) |
| Pre-defined operating mode | YES (supervised by the controller according to the user-defined ozone concentration and sterilization time) |
| Self-diagnosis | YES |
| Working time counter | YES (signaling the replacement of filters and UV radiators) |
| Power supply | 230 V (AC), 50 Hz |
| Rated current | 2.0 A |
| Rated power | 470 W |
| Power cable length | 3 m / 10 m* |
| Dimensions (H x W x D) | 1490 x 700 x 760 |
| Net weight | 26 kg |
| Type of housing | metal, powder-coated |
| Transport wheels/handles | transport wheels |
| Additional functions | signaling safe and excessive ozone levels, detachable power cord |

Recommended
CUBIC CAPACITY
of the sterilized room

850 m³

WORK MODES:



Filtration



Disinfection
UV-C



UV-C
disinfection -
silent mode




Sterilization
O₃ + UV-C (Auto)



TECHNICAL SPECIFICATIONS:

| MODE – STERILIZATION O ₃ + UV-C | |
|---|-----------------------|
| Recommended maximum cubic capacity of the sterilized room (to achieve 5 ppm O ₃) | 850 m ³ |
| Ozone generator capacity | 30,000 mg/h |
| Fan capacity Sterilization mode O ₃ +UV-C | 410 m ³ /h |
| Noise level – sterilization mode O ₃ +UV-C | 55 dB(A) |
| STANDARD MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 575 J/m ² |
| Average UV-C radiation dose | 141 J/m ² |
| Fan performance UV-C disinfection mode | 290 m ³ /h |
| Noise level – disinfection mode UV-C | 49 dB(A) |
| SILENT MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 901 J/m ² |
| Average UV-C radiation dose | 221 J/m ² |
| Fan performance UV-C disinfection mode (silent) | 185 m ³ /h |
| Noise level – disinfection mode UV-C (silent) | 33 dB(A) |
| FILTRATION MODE | |
| Fan capacity in filter mode | 290 m ³ /h |
| Noise level – filtering mode | 49 dB(A) |
| UV lamp | YES |



| | |
|------------------------------|---|
| Type lamps UV | UV-C sterilization λ = 253.7 nm |
| Electrical power of UV lamps | 330 W |
| Durability of UV lamps | 9,000 h |
| Power of the UV light source | 118 W |
| Ozone generator | YES |
| Ozone concentration sensor | YES |
| Ozone destructor | YES |
| Air filtering | 2-stage |
| Fan type | Radial |
| Fan engine | enclosed |
| Type | flowtype |
| Controller | YES |
| Automatic mode of operation | YES (control of sterilization time according to actual ozone concentration readings) |
| Pre-defined operating mode | YES (supervised by the controller according to the user-defined ozone concentration and sterilization time) |
| Self-diagnosis | YES |
| Working time counter | YES (signaling the replacement of filters and UV radiators) |
| Power supply | 230 V (AC), 50 Hz |
| Rated current | 2.9 A |
| Rated power | 670 W |
| Power cable length | 3 m / 10 m* |
| Dimensions (H x W x D) | 1490 x 700 x 760 |
| Net weight | 36 kg |
| Type of housing | metal, powder-coated |
| Transport wheels/handles | transport wheels |
| Additional functions | signaling safe and excessive ozone levels, detachable power cord |

Recommended
CUBIC CAPACITY
of the sterilized room

850 m³

WORK MODES:



Filtration



Disinfection
UV-C



UV-C
disinfection -
silent mode





Sterilization
O₃ + UV-C (Auto)



TECHNICAL SPECIFICATIONS:

| | |
|---|-----------------------|
| MODE – STERILIZATION O₃ + UV-C | |
| Recommended maximum cubic capacity of the sterilized room (to achieve 5 ppm O ₃) | 850 m ³ |
| Ozone generator capacity | 30,000 mg/h |
| Fan capacity Sterilization mode O ₃ +UV-C | 900 m ³ /h |
| Noise level – sterilization mode O ₃ +UV-C | 61 dB(A) |
| STANDARD MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 359 J/m ² |
| Average UV-C radiation dose | 90 J/m ² |
| Fan performance UV-C disinfection mode | 680 m ³ /h |
| Noise level – disinfection mode UV-C | 58 dB(A) |
| SILENT MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 421 J/m ² |
| Average UV-C radiation dose | 105 J/m ² |
| Fan performance UV-C disinfection mode (silent) | 580 m ³ /h |
| Noise level – disinfection mode UV-C (silent) | 55 dB(A) |
| FILTRATION MODE | |
| Fan capacity in filter mode | 680 m ³ /h |
| Noise level – filtering mode | 58 dB(A) |
| UV lamp | |
| YES | |



| | |
|------------------------------|---|
| Type lamps UV | UV-C sterilization λ = 253.7 nm |
| Electrical power of UV lamps | 440 W |
| Durability of UV lamps | 9,000 h |
| Power of the UV light source | 157 W |
| Ozone generator | YES |
| Ozone concentration sensor | YES |
| Ozone destructor | YES |
| Air filtering | 2-stage |
| Fan type | Radial |
| Fan engine | enclosed |
| Type | flowtype |
| Controller | YES |
| Automatic mode of operation | YES (control of sterilization time according to actual ozone concentration readings) |
| Pre-defined operating mode | YES (supervised by the controller according to the user-defined ozone concentration and sterilization time) |
| Self-diagnosis | YES |
| Working time counter | YES (signaling the replacement of filters and UV radiators) |
| Power supply | 230 V (AC), 50 Hz |
| Rated current | 3.5 A |
| Rated power | 810 W |
| Power cable length | 3 m / 10 m* |
| Dimensions (H x W x D) | 1490 x 700 x 760 |
| Net weight | 48 kg |
| Type of housing | metal, powder-coated |
| Transport wheels/handles | transport wheels |
| Additional functions | signaling safe and excessive ozone levels, detachable power cord |

Recommended
CUBIC CAPACITY
of the sterilized room

1050 m³

WORK MODES:



Filtration



Disinfection
UV-C



UV-C
disinfection -
silent mode




Sterilization
O₃ + UV-C (Auto)



TECHNICAL SPECIFICATIONS:


| | |
|--|-------------------------|
| MODE – STERILIZATION O₃ + UV-C | |
| Recommended maximum cubic capacity of the sterilized room (to achieve 5 ppm O ₃) | 1,050 m ³ |
| Ozone generator capacity | 40,000 mg/h |
| Fan capacity Sterilization mode O ₃ +UV-C | 1,180 m ³ /h |
| Noise level – sterilization mode O ₃ +UV-C | 57 dB(A) |
| STANDARD MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 266 J/m ² |
| Average UV-C radiation dose | 79 J/m ² |
| Fan performance UV-C disinfection mode | 900 m ³ /h |
| Noise level – disinfection mode UV-C | 54 dB(A) |
| SILENT MODE – UV-C DISINFECTION | |
| Initial maximum UV-C radiation dose | 352 J/m ² |
| Average UV-C radiation dose | 104 J/m ² |
| Fan performance UV-C disinfection mode (silent) | 680 m ³ /h |
| Noise level – disinfection mode UV-C (silent) | 52 dB(A) |
| FILTRATION MODE | |
| Fan capacity in filter mode | 900 m ³ /h |
| Noise level – filtering mode | 54 dB(A) |
| UV lamp | YES |



| | |
|------------------------------|---|
| Type lamps UV | UV-C sterilization λ = 253.7 nm |
| Electrical power of UV lamps | 550 W |
| Durability of UV lamps | 9,000 h |
| Power of the UV light source | 196 W |
| Ozone generator | YES |
| Ozone concentration sensor | YES |
| Ozone destructor | YES |
| Air filtering | 2-stage |
| Fan type | Radial |
| Fan engine | enclosed |
| Type | flowtype |
| Controller | YES |
| Automatic mode of operation | YES (control of sterilization time according to actual ozone concentration readings) |
| Pre-defined operating mode | YES (supervised by the controller according to the user-defined ozone concentration and sterilization time) |
| Self-diagnosis | YES |
| Working time counter | YES (signaling the replacement of filters and UV radiators) |
| Power supply | 230 V (AC), 50 Hz |
| Rated current | 4.4 A |
| Rated power | 1,010 W |
| Power cable length | 3 m / 10 m* |
| Dimensions (H x W x D) | 1490 x 700 x 760 |
| Net weight | 55 kg |
| Type of housing | metal, powder-coated |
| Transport wheels/handles | transport wheels |
| Additional functions | signaling safe and excessive ozone levels, detachable power cord |

Bibliography:

| | |
|------|--|
| [1] | UV doses from the publication "Disinfection by UV-radiation" by PHILIPS |
| [2] | "UVC LED Irradiation Effectively Inactivates Aerosolized Viruses, Bacteria, and Fungi in a Chamber-Type Air Disinfection System"; Do-Kyun Kim, Dong-Hyun Kang; August 2018; American Society for Microbiology Journals |
| [3] | "2020 COVID-19 Coronavirus Ultraviolet Susceptibility"; W. J. Kowalski, T.J Walsh, V. Petraitis, March 2020, ResearchGate |
| [4] | www.clordisys.com/pdfs/misc/UV%20Data%20Sheet.pdf |
| [5] | www.boviemedical.com/wp-content/uploads/2018/04/uv24-lab-results-kowalski-wp-aerobiology.pdf |
| [6] | "Ozonoterapia oraz zastosowanie ozonu w dezynfekcji"; D. Białoszewski, E. Bocian, S. Tyski, May 2020, POST. MIKROBIOL. 2012, 51, 3, 177–184 |
| [7] | "Use of ozone in the food industry"; Zeynep B. Guzel-Seydim, Annel K. Greene, A.C. Seydim, LWT - Food Science and Technology, Volume 37, Issue 4, June 2004, Pages 453–460 |
| [8] | "Ozonation and UV irradiation – an introduction and examples of current applications"; Steven T. Summerfelt, Aquacultural Engineering, Volume 28, Issues 1–2, June 2003, Pages 21–36 |
| [9] | "Zastosowanie ozonu w przemyśle spożywczym"; K. Krosowiak, K. Śmigielski, P. Dziugan, Przemysł Spożywczy 11/2017 |
| [10] | "Fluence (UV Dose) Required to Achieve Incremental Log Inactivation of Bacteria, Protozoa, Viruses and Algae"; Adel Haji Malayeri, Madjid Mohseni, Bill Cairns, James R. Bolton, Gabriel Chevretils, Eric Caron, 2006 |
| [11] | "Ozone Disinfection of SARS–Contaminated Areas"; Kenneth K. K. LAM |
| [12] | "Ozonation and UV Disinfection"; Steven Summerfelt & Brian Vinci; Freshwater Institute, Shepherdstown, WV |
| [13] | "Molecular Mechanisms of Ultraviolet Radiation–Induced DNA Damage and Repair"; R. P. Rastogi, Richa, A. Kumar, M.B. Tyagi, R.P. Sinha; Journal of Nucleic Acids, Volume 2010 |
| [14] | "UVC photon-induced denaturing of DNA: A possible dissipative route to Archean enzyme-less replication"; Karo Michaelian, Norberto Santillan Padilla; Heliyon, Volume 5, Issue 6, June 2019, e019025 |



PROFESSIONAL ROOM STERILIZERS

STERYLIS®

ULTRA



MILOO-ELECTRONICS Sp. z o. o.
Stary Wiśnicz 289
32-720 Nowy Wiśnicz
Poland



www.sterylis.com